


AlloDerm® Regenerative Tissue Matrix



**Soft tissue replacement
without a palatal harvest**

BIOHORIZONS®
SCIENCE • INNOVATION • SERVICE



AlloDerm



What is AlloDerm Regenerative Tissue Matrix?

AlloDerm is an acellular dermal matrix derived from donated human skin that undergoes a multi-step proprietary process that removes both the epidermis and the cells that can lead to tissue rejection. AlloDerm has been used in a wide variety of soft tissue grafting procedures such as root coverage, soft tissue augmentation and guided bone regeneration with a consistent record of excellent results.¹⁻⁷



AlloDerm offers numerous advantages compared to the connective tissue autograft from the patient's palate:

- Eliminates the need for palatal surgery
- Removes palatal harvesting limitations from treatment planning considerations
- Reduces patient reluctance to follow through with surgical treatment
- Consistent quality
- Provided in multiple convenient sizes
- Available in two thickness ranges for use in different procedures:
 - 0.9 to 1.6 mm - AlloDerm for root coverage, soft tissue ridge augmentation, etc.
 - 0.5 to 0.8 mm - AlloDerm GBR for guided bone regeneration and barrier membrane function

Procurement and safety

AlloDerm has a safety history of more than a decade. Introduced in 1994 for treating burn patients, AlloDerm has proved its versatility and safety in more than a million diverse procedures in general, orthopedic, urogenital, and dental surgeries.⁸

AlloDerm owes its exemplary safety to the safeguards at every step starting from donor screening to the final packaging.

- Tissue accepted only from AATB (American Association of Tissue Banks) compliant tissue banks
- Extensive panel of serology tests
- Proprietary processing technology removes immunogenic cells and minimizes risk of disease transmission
- Final sterility testing ensures that no external pathogens are introduced while processing

Regenerative Tissue Matrix


Processing of AlloDerm

The proprietary processing to derive AlloDerm from donor tissue involves a series of steps:

- Treatment with buffered salt solution to separate and eliminate the epidermis
- Series of washes with mild non-denaturing detergent solutions to solubilize and eliminate all cells
- Final freeze drying step using patented technology that prevents damaging ice crystal formation


Allograft Tissue

Removal of epidermis and cells



Regenerative Tissue Matrix

Complex acellular heterogenous scaffold, with growth factor binding sites and blood vessel architecture; dehydrated and ready to implant




How does AlloDerm work?

AlloDerm provides a matrix consisting of collagens, elastin, vascular channels, and proteins that support revascularization, cell repopulation and tissue remodeling.


After placement, the patient's blood infiltrates the AlloDerm graft through retained vascular channels, bringing host cells that adhere to proteins in the matrix. Significant revascularization can begin as early as one week after implantation. The host cells respond to the local environment and the matrix is remodeled into the patient's own tissue, in a fashion similar to the body's natural tissue attrition and replacement process.



Because the components remain in their natural biologically active state, ADM is immediately recognized as human tissue.



Rapid cell repopulation and revascularization. Initiation of intrinsic regeneration process.



Complete remodeling into the patient's own tissue. Functional, physiological and reconstructive outcomes.

AlloDerm

Documented equivalence to autogenous connective tissue


Multiple, randomized clinical trials (RCT) have shown root coverage results with AlloDerm to be equivalent to autogenous connective tissue, and concluded that the procedure was predictable and practical. A meta-analysis of eight RCTs showed no statistically significant differences between the two groups for measured outcomes: recession coverage, keratinized tissue formation, probing depth and clinical attachment levels.⁹

Keys for successful root coverage include:


- Thorough root conditioning and/or restoration removal
- Flap or Pouch design that minimizes loss of vascularity
- Tension-free coronal positioning of flap or pouch to completely cover AlloDerm

Histological evidence of remodelling

A human histologic evaluation of AlloDerm and connective tissue (CT) documented that both formed a band of dense collagenous tissue when placed beneath a coronally advanced flap. Gingival attachment, a combination of long junctional epithelium and connective tissue adhesion, was comparable for both groups. At six months postoperatively, the overall histologic outcomes were similar for both CT and AlloDerm grafts.¹⁰



Connective tissue specimen demonstrating mucosal tissue (M) overlying dense grafted connective tissue (C) and osseous crest (B). Original magnification 40X; hematoxylin and eosin (H&E) staining.



Acellular dermal matrix specimen demonstrating mucosal tissue (M) overlying the area of graft placement (ADM) and osseous crest (B). Original magnification 40X. Verhoff solution stained elastin fibers help differentiate graft area from host tissue.

Cummings et al. Histologic Evaluation of Autogenous Connective Tissue and Acellular Dermal Matrix Grafts in Humans. J Periodontol 2005;76:178-186.

Ease of use




Remove from pouch

- AlloDerm has up to a 2-year shelf-life when stored between 1°-10°C (34°-50°F).
- Open outer foil pack. Drop graft into saline bath directly from inner pouch.



Rehydrate

- Rehydrate in two consecutive sterile saline baths.
- Remove paper backing from AlloDerm between first and second baths.



Two distinct surfaces

- Basement membrane side (BM) is rough and will not readily absorb blood.
- Dermal side is smooth and will absorb blood.

Important: Before use, clinicians should review all risk information and directions, which can be found on the packaging and in the "Information for Use" attached to the packaging of each AlloDerm graft.

Regenerative Tissue Matrix

Root Coverage

AlloDerm is ideal for treating multiple defects in a single procedure. Available sizes include: 1cm x 1cm, 1cm x 2cm, 1cm x 4cm and 2cm x 4cm. After hydration it may be trimmed to the desired size with a scalpel or sharp scissors.

Photos courtesy of Dr. Edward P. Allen, Dallas, Texas



Gingival recession with root surface restorations.



AlloDerm graft placed in pouch and sutured.



Complete root coverage at one year postoperatively.

Soft Tissue Ridge Augmentation

AlloDerm can be used effectively for soft tissue ridge augmentation. A tunnel or pouch may be created beneath the defect into which the AlloDerm can be inserted. If multiple layers of AlloDerm are used for increased thickness, it is recommended that it be layered, rather than rolled. In this indication, orient the dermal surfaces on the outside of the graft.

Photos courtesy of Dr. Edward P. Allen, Dallas, Texas



Alveolar ridge deficiency at site of missing maxillary left lateral incisor and canine.



Folded and sutured AlloDerm graft placed and sutured within the soft tissue pouch.



3 months post-op showing restoration of normal alveolar ridge contour.

Soft Tissue Augmentation Around Dental Implants

AlloDerm is effective in augmenting thin tissue around dental implants to create more attached tissue.

Photos courtesy of Dr. Carl E. Misch, Beverly Hills, Michigan



Treatment plan for revision of a failing 2-implant overdenture to a 5-unit cemented bar overdenture.



The bony defect is grafted using autologous bone from osteotomies. The AlloDerm is oriented with basement membrane up.



Postoperative results show thick, immobile tissue.

A biological barrier for better hard and soft tissue regeneration



Unlike other barrier membranes that either resorb too quickly or do not resorb at all, AlloDerm GBR actually allows the body to remodel it into the patient's own tissue. This results in good bone regeneration while benefiting soft tissue quality and esthetics.¹¹⁻¹²

AlloDerm GBR, unlike conventional barrier membranes, if left exposed because of incomplete primary closure, has a significantly reduced chance of infection and graft failure due to rapid revascularization. A recent study demonstrated that sites covered with AlloDerm generated 16% more vital bone in extraction sockets than did sites covered with an ePTFE membrane.¹³


A separate study concluded that AlloDerm, used as a barrier over resorbable hydroxyapatite in extraction sites, was able to preserve ridge dimensions and significantly increase the width of keratinized tissue.¹⁴

Guided Bone Regeneration

AlloDerm GBR is produced by exactly the same process as AlloDerm, but it is thinner (0.5 to 0.9mm) which facilitates primary closure. It readily adapts to graft sites and can be secured with either sutures or tacks.



Extraction sites grafted with GRAFTON® DBM.



AlloDerm GBR covering the ridge beneath the flap to obtain tension-free closure.




Post-operative result shows excellent soft tissue coverage.

Sizes available




2cm x 4cm



2cm x 2cm



1cm x 4cm



1cm x 2cm



1cm x 1cm

AlloDerm and AlloDerm GBR are available in these convenient sizes. Examples are shown to scale.

Regenerative Tissue Matrix

References

1. Management of Gingival Recession by the Use of a Acellular Dermal Graft Material: A 12-Case Series. Santos A, Goumenos G and Pascual A. *J Periodontol* 2005;76(11):1982-1990.
2. Subpedicle Acellular Dermal Matrix Graft and Autogenous Connective Tissue Graft in the Treatment of Gingival Recessions: A Comparative 1-Year Clinical Study. Paolantonio M, Dolci M, Esposito P, D'Archivio D, Lisanti L, Di Luccio A and Perinetti G. *J Periodontol* 2002;73(11):1299-1307.
3. Clinical Evaluation of Acellular Allograft Dermis for the Treatment of Human Gingival Recession. Aichelmann-Reidy ME, Yukna RA, Evans GH, Nasr HF and Mayer ET. *J Periodontol* 2001;72(8):998-1005.
4. Predictable Multiple Site Root Coverage Using an Acellular Dermal Matrix Allograft. Henderson RD, Greenwell H, Drisko C, Regennitter FJ, Lamb JW, Mehlbauer MJ, Goldsmith LJ and Rebitski G. *J Periodontol* 2001;72(5):571-582.
5. Surgical therapies for the treatment of gingival recession. A systematic review. Oates TW, Robinson M and Gunsolley JC. *Ann Periodontol* 2003;8:303-320.
6. Root coverage of advanced gingival recession: A comparative study between acellular dermal matrix allograft and subepithelial connective tissue grafts. Tal H, Moses O, Zohar R, et al. *J Periodontol* 2002;73:1405-1411.
7. The clinical effect of acellular dermal matrix on gingival thickness and root coverage compared to coronally positioned flap alone. Woodyard JG, Greenwell H, Hill M, et al. *J Periodontol* 2004;75:44-56.
8. Data on file.
9. Acellular Dermal Matrix for Mucogingival Surgery: A Meta-Analysis. Gapski R, Parks CA and Wang HL. *J Periodontol* 2005;76(11):1814-1822.
10. Histologic Evaluation of Autogenous Connective Tissue and Acellular Dermal Matrix Grafts in Humans. Cummings LC, Kaldahl WB and Allen EP. *J Periodontol* 2005;76(2):178-186.
11. Acellular Dermal Matrix Graft as a Membrane for Guided Bone Regeneration: A Case Report. Novaes AB Jr and Souza SL. *Implant Dentistry* 2001;10(3):192-195.
12. Ridge Preservation Utilizing an Acellular Dermal Allograft and Demineralized Freeze-Dried Bone Allograft: Part I. A Report of 2 Cases. Fowler EB, Breault LG and Rebitski G. *J Periodontol* 2000;71(8):1353-1359.
13. Extraction Sockets and Implantation of Hydroxyapatites With Membrane Barriers, A Histologic Study. Froum S, Cho SC, Elian N, Rosenberg E, Rohrer M and Tarnow D. *Implant Dentistry* 2004;13(2):153-164.
14. Acellular Dermal Matrix and Hydroxyapatite in Prevention of Ridge Deformities after Tooth Extraction. Luczyszyn SM, Papalexiou V, Novaes AB Jr, Grisi MF, Souza SL and Taba M Jr. *Implant Dent* 2005;14(2):176-184.

GRAFTON® DBM and MINERoss® are processed by:



AlloDerm® and AlloDerm® GBR™ are processed by:



BioHorizons USA

2300 Riverchase Center
Birmingham, AL 35244

Customer Care / Servicio al Cliente:
888-246-8338 or 205-967-7880

BioHorizons Canada

21 Amber Street, Unit # 6
Markham, Ontario L3R 4Z3

Customer Care / Service à la Clientèle:
866-468-8338 or / ou 905-944-1700

BioHorizons Spain

Serrano Anguita, 10
28004 Madrid, España

Atención al Cliente:
+34 91 713 10 84

BioHorizons UK

17 Wellington Business Park
Dukes Ride
Crowthorne, Berkshire RG45 6LS

Customer Care:
+44 1344 752560

BioHorizons Germany

Marktplatz 3
79199 Kirchzarten

Kunden Service:
+49 7661-909989-0

BioHorizons Australia

25-33 Allen Street
Waterloo NSW 2012

Customer Care:
+61 2 8399 1520

BioHorizons Mexico

Kelvin 8 Dept. 303
Col. Anzures
C.P. 11590, Mexico, D.F.

Servicio al Cliente:
+52 55 5545 1297

BioHorizons Chile

Cerro Colorado 5030, Oficina 513
Las Condes
Santiago, Chile

Atención al Cliente:
+56 2 361 9519

As applicable, BioHorizons products are cleared for sale in the European Union under the EU Medical Device Directive 93/42/EEC. We are proud to be registered to ISO 13485:2003, the international quality management system standard for medical devices, which supports and maintains our product licences with Health Canada and in other markets around the globe.

Not all products are available in all markets. AlloDerm® and AlloDerm® GBR™ must be shipped overnight and may not be returned or exchanged for credit due to specific guidelines regarding the storage of allograft tissue.

© 2008 BioHorizons Implant Systems, Inc. All Rights Reserved. MLD102 REV A NOV 2008

w w w . b i o h o r i z o n s . c o m

